FMA - Functional Mandibular Advancer*

Ref.-Nr. 330-0100

GB *Laboratory instructions and
Doctor's instructions*

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CE 0297 93/42/EWG

FORESTADENT®
GERMAN PRECISION IN ORTHODONTICS

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Disclaimer

The Forestadent Bernhard Förster GmbH product guarantee covers material or manufacturing faults only. The operator bears sole responsibility for correct handling and use of products. The Forestadent guarantee does not cover defects that occur due to inadequate treatment planning, incorrect fabrication or use of the treatment appliance.

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
Option C • Variante C

GB

- a) Adapt bands intraorally and transfer them by means of an alginate impression to the working model. Bands may be used on first premolar and first molar (option A) as well as on first and second molar (option B).
- b) Adapt transpalatal and lower lingual archwires or use a Slim Line expander for expansion.
- c) If the second molars are only partially erupted, occlusal rests can be placed to ensure stress distribution (Option C).

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GB


- a) Take a facebow registration and mount the upper jaw model into the articulator.
- b) Adjust lower jaw model 3-dimensionally into the therapeutic position with help of the construction bite. The more distally the FMA is fitted (first and second molars), the greater the intraoral comfort. The FMA should be placed as far dorsally as possible. Brackets can be fitted to the dental arch from teeth 5 – 5 to perform simultaneous treatment options. If the second molar has not yet erupted, the first premolar should be incorporated. Ensure adequate occlusal support!
- c) Hold the protrusive bars ① and the inclined planes ③ in contact and adjust them parallel to the occlusal plane. The protrusive bars must be aligned according the articulation excursions in a programmed articulator to improve the intraoral comfort.

- d) Fix the upper and lower mounting plates ② & ③ to the bands with sticky wax or temporary spot welding.
- e) Secure the protrusive bars ① into the middle position on the mounting plate to gain more flexibility when activating (see 2).

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3 Laboratory instructions • Laboranleitung

GB



- a) Cut down the mounting plates ② & ③ if required. Check the position of all FMA parts in the articulator in relation to the joints by simulating the masticatory movements.
- b) Simulate all lower jaw movements and adjust the protrusive bars ① in relation to the lower vestibulum.
- c) Weld the mounting plates ② & ③ onto the bands in their final position.

The planned occlusion must be easily repositionable on the stone model. Open the bite approx. 1 mm in the articulator! The physiological temporomandibular joint relationships of the planned occlusion should be checked in an MRT, at least in the case of adults. Metal spurs must be interference free during protrusion and laterotrusion in opening and closing movements.

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GB

- a) Tack weld the palatal and the lingual archwire to the bands.
- b) Apply flux on mounting plates ② & ③, bands, transpalatal and lower lingual archwire and weld all parts very carefully to avoid repairs (parts can also be lasered). Metal joints are temperature resistant up to approx. 700°C; the bars should therefore not be cast on. The existing solder joints between the bar and body must be covered with heat barrier paste. Cast bands should be as wide as possible, to enlarge the adhesive surface and increase stability.

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GB

Polish and check all parts, mount the appliance in the articulator and recheck all lower jaw movements in the therapeutic position and cut down the protrusive bars ① if necessary. The use of dual-curing cement and conditioning of the tooth improves the bond strength.


A guidance extension can be welded onto the FMA, if there is a risk of disengaging due to lateral clearance. ①

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**GB**

Ensure that the retention screws are fixed securely in position! The screws can be fixed in position using standard techniques, e.g. by applying composite.

It is important to follow the stages in the instructions for use exactly when fabricating and using the appliance. Courses offered by Forestadent provide support in the correct fabrication of the appliance. We recommend attending these courses.

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GB**After fitting:**

- Protrusive/ opening exercises should be carried out with the patient.
- The patient should be fully informed about the risks with fixed appliances.
- The patient should be instructed on the necessity of avoiding sticky/ hard foods.

Allow an acclimatisation period of one week. Check for correct adaptation no later than two weeks after fitting.

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