Frog II Appliance Laboratory instructions*
Frosch II Apparatur Laboranleitung*

Ord.-Nr. A167D1639

CE 0297

FORESTADENT®
GERMAN PRECISION IN ORTHODONTICS

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
The Frog II screw has two long retention arms on one side. These can be bent and therefore easily adjusted to the position of the abutments of the OrthoEasy® Pal system.

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Upper Fig.: Impression with embedded impression caps (Order-no.1299-1005). After cleaning the impression check if the impression caps are correctly embedded. Place the lab analogs (Order No. 1299-1001) in the impression caps.

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
The lab analogs can be slightly waxed in. Then cast the plaster model with super hard stone (class IV).

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
Position of the Ortho Easy Pal pin is duplicated with the lab analogs in the plaster model. The OrthoEasy Pal abutments (Order No.1299-1004) are placed on the lab analogs. To keep them in the correct position they should be fixated with the retaining screws for the palatal pin (order-no. 1299-1003). This prevents the abutments from moving when the wire is bent around them (see next step).

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Once the abutments are placed on the lab analogs and fixed with retaining screws, bend the retention arms of the Frog II screw around them. This should be done in a way so that the Frog II screw can be correctly placed in the palate. Always use two pliers for bending. It should be noted that the weld seam is protected. It is recommended to bend the retention arms with two pliers. When bending the retention arms please make sure that you do not put any pressure on the welded joint.

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The screw is placed on the level of the center of resistance of the molars. Following points have to be considered:

- Retention arms must not touch the soft tissue.
- The Frog II screw has to be placed parallel to the occlusion plane and as close to the level of the center of resistance of the molars as possible.
- Fixation of the screw on the model is done with silicone.

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After the silicone has set the retention arms are bent around the abutments and shortened appropriately. Now the retention arms can be spot or laser-welded. Laser-welding can be done on the model. With spot-welding it is recommended to solder subsequently. To do this, take the Frog II screw with the spot-welded abutments off the model. Make sure that no solder flows into the abutments.

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
Wire bending

Use only TMA wires with \varnothing 0,8 mm. Start with the loop which will be inserted into the lock on the lingual side. Now make a U-bend (1) in dorsal direction with a double-bend (2) at the end. Repeat the same on the other side. Make sure that the double-bend is placed slightly more distal than the mesial opening of the lingual sheath on the molar band.

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As shown the double-bend is placed distally of the mesial opening of the lingual sheath (1). Additionally it is bend by around 15° towards occlusal. Check the force of the spring with a stress and tension gauge. The force exerted by the spring should be at least 180 cN (2).

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Following points have to be considered during spring fabrication:

- U-bend points to dorsal
- The double-bend is placed slightly distal to the mesial opening of the lingual sheath of the molar band
- To avoid distal tipping of the molar the double-bend is bent by around 15° towards occlusal
- Attention: The double-bend must be level with the band sheath. A position too low leads to intrusion. A position too high leads to extrusion.
- When inserting the double-bends into the lingual sheath the spring should exert a force of at least 180 cN. This force shows how much the spring is activated distally. Should it be too low the double-bend has to be bent more towards dorsal or the U-bends have to be made smaller.
- The fabricated spring must not touch any soft tissue when inserted and at any time during treatment.
- Finally chamfer the ends of the double-bends and bend them towards the teeth as much as possible to avoid injuries to the tongue.

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Phone +49 (0) 7231 459-0 · Fax +49 (0) 7231 459-102

info@forestadent.com · www.forestadent.com

FORESTADENT[®]
GERMAN PRECISION IN ORTHODONTICS

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Bernhard Förster GmbH